

# MarvelXACT™ PEEKsil™

## Product Data Sheet



Fused  
Silica Tubing



Ultra-high  
Pressure



Torque Limiting  
Technology



Learn More at  
[idex-hs.com](http://idex-hs.com)



100%  
Bioinert



Our bioinert MarvelXACT PEEKsil connection systems utilize patented MarvelXACT fittings and PEEKsil tubing for use in the most demanding UHPLC applications. The MarvelXACT fitting features torque-limiting technology that is finger tight and provides a haptic “click” feedback when it reaches the optimum torque, assuring a perfect zero dead volume connection every time. The robust PEEKsil tubing has smooth silica glass wetted surfaces and is chemically compatible to organic solvents. Each connection system is individually cut, polished, and manufactured to provide industry leading technology that produces superior chromatographic results.

FINGER-TIGHT TO  
19,000 psi

REUSABLE  
100+  
USES

FLAT-FACED  
SEALING CONNECTIONS  
ZERO  
DEAD VOLUME

