


Nanopede™ Filter Set

FLIGHT 2

12.5 mm




For Technical Questions:
Semrock@idexcorp.com

Semrock

FLIGHT 2 12.5 mm

■ UV ■ VIS ■ NIR



Scan to Visit our eStore
for more Nanopede
options, Technical
Information, and
Specifications:

index-hs.com/flight2-12_5

C-mount

C-mount

C-mount

FBP01-350/17-12.5

FBP01-370/16-12.5

FBP01-390/16-12.5

FBP01-410/16-12.5

FBP01-430/16-12.5

FBP01-450/16-12.5

FBP01-470/15-12.5

FBP01-490/15-12.5

FBP01-510/15-12.5

FBP01-530/15-12.5

FBP01-550/15-12.5

FBP01-570/14-12.5

FBP01-590/14-12.5

FBP01-610/14-12.5

FBP01-630/14-12.5

FBP01-650/14-12.5

FBP01-670/13-12.5

FBP01-690/13-12.5

FBP01-710/13-12.5

FBP01-730/13-12.5

FBP01-750/13-12.5

FBP01-775/22-12.5

FBP01-805/22-12.5

FBP01-835/22-12.5

FBP01-865/21-12.5

FBP01-895/21-12.5